







# Orthomed

Manufacturer and Supplier of  
High Quality, Innovative Veterinary  
Orthopaedic Products

# About Orthomed

## Meeting and exceeding medical gold standards

For over ten years, we have been developing, manufacturing and distributing a range of technically advanced implants and instruments which are engineered to simplify the most complex surgeries, underpinned by service, support and training which empowers specialists, general practitioners and first opinion vets to carry out gold standard surgical procedures.

Our product portfolio includes some of the most innovative orthopaedic products and systems available to the veterinary industry: including products such as MMP OrthoFoam™, SOP™ and RidgeStop™. Product development is carried out by bio-mechanical engineers and material scientists, who work with world renowned orthopaedic surgeons, all driven by the desire to bring to market products which break new ground whilst meeting and exceeding medical gold standards for surgical equipment.

We believe that the sale is just the first part of the journey; our customers are supported for life by a multi-faceted service offering, which includes hands-on surgical training workshops, pre and post-surgery advice and support, and access to a wide range of educational tools.

We are inspired by ground breaking engineering and satisfied only by enabling the best possible outcome for surgeon, patient and owner. Orthomed is a truly global company, headquartered in the UK, with offices in Europe, Australasia, North America, South Africa and the Far East.

- ✓ **Leaders in Innovation**
- ✓ **Surgical Training Workshops**
- ✓ **Aftersales Support**
- ✓ **High Quality Products**
- ✓ **Backed by Industry Experts**
- ✓ **Advanced Ongoing R&D**



# Contents

---

<b>Workshops</b>	<b>6</b>
01 - <b>MMP</b>	<b>18</b>
02 - <b>RidgeStop™</b>	<b>28</b>
03 - <b>SOP™</b>	<b>36</b>
04 - <b>Carpal Arthrodesis</b>	<b>50</b>
05 - <b>TPLO</b>	<b>54</b>
06 - <b>Internal Fixation</b>	<b>58</b>
07 - <b>Instruments</b>	<b>72</b>
08 - <b>Power Tools</b>	<b>80</b>
<b>More From Orthomed</b>	<b>84</b>
<b>Orthomed around the world</b>	<b>88</b>







# Workshops

We are not just about selling products

---






Our aim is to inspire specialists, first opinion vets and general practitioners to be able to carry out gold standard surgical procedures, and we aim to do this through education and teaching.




All of our workshops consist of a mixture of theory, followed by a practical session using implants on sawbones which can then be taken away to show prospective patient owners. We also run cadaver workshops throughout the year.



We do this to ensure that the surgeon is 100% certain that they are confident in the system, that they understand the technique, the planning, the approach, the aftercare as well as avoiding any potential pitfalls. We want all surgeons to go away knowing that they have made the right choice.



### Order by Phone

 +44 (0) 845 045 0259  
 +1 772-562-6044  
 +61 (0) 8 9590 8850

 +49 (0) 4191 8030013  
 +34 902 47 47 01  
 +39 (0) 391 11 60 784

 +27 (0) 83 227 8181  
 +886 4 2258 6668

### Order Online





All of our products are available to buy online at:

[www.orthomed.co.uk/nz](http://www.orthomed.co.uk/nz)

# Workshops


Global workshops held all year round...




-  **Presented by industry experts and diplomats**
-  **Wet labs and saw-bone workshops**
-  **Held regularly world-wide**
-  **Gives you complete confidence to perform the surgery**







**Order by Phone**
 +44 (0) 845 045 0259

 +1 772-562-6044


 +61 (0) 8 9590 8850

 +49 (0) 4191 8030013

 +34 902 47 47 01

 +39 (0) 391 11 60 784

 +27 (0) 83 227 8181

 +886 4 2258 6668
**Order Online**

All of our products are available to buy online at:

[www.orthomed.co.uk/nz](http://www.orthomed.co.uk/nz)

9

**About our Workshops**


We are aware that the success of our products hinges upon the experiences of the surgeons who use them, and with this in mind – educating them in the use of the products is a hugely important part of what we do.

We have worked with leading experts in the field of veterinary orthopaedics on the development of the products, and the knowledge and experience gained from these partnerships is passed on via our hands-on workshops.

Workshops take place throughout the year, across the globe, and are pitched at levels suitable for surgeons at every stage of their professional development. They are classroom based and offer participants the chance to carry out procedures, such as RidgeStop or MMP, on plastic models or cadavers, under the instruction of highly regarded tutors.

**Our Workshops...**


Details of our courses are available on the following pages or on the Orthomed website. [www.orthomed.co.uk/nz](http://www.orthomed.co.uk/nz)




MMP



RidgeStop™



SOP™



Advanced Spinal Course

Workshops

# Course Outlines


“


*Great and prompt service. Very knowledgeable, effective service, next day delivery and high quality implants. Loving the MMP technique!*”

**Eric Rouviere**

Rotherham Vets4Pets


---

**Order by Phone**
 +44 (0) 845 045 0259

 +1 772-562-6044


 +61 (0) 8 9590 8850

 +49 (0) 4191 8030013

 +34 902 47 47 01

 +39 (0) 391 11 60 784

 +27 (0) 83 227 8181

 +886 4 2258 6668
**Order Online**

All of our products are available to buy online at:

[www.orthomed.co.uk/nz](http://www.orthomed.co.uk/nz)

11

# MMP

## TTA for Canine Cruciate Repair

### Course Content

The Modified Maquet Procedure (MMP) is an evolution of the TTA procedure for treatment of lameness due to cranial cruciate failure. The course is in two parts: a lecture based introduction followed by a practical session using synthetic bone models of the stifle, animations and video support material will also be used. The development of the technique, implants and associated instruments will be described and a detailed review of surgical technique and clinical outcomes will cover all aspects of MMP.


### Benefits

- Over 35,000 dogs treated\*
- Over 800 users globally\*
- Use of a load-bearing, bio-compatible pure titanium wedge
- Rapid bone ingrowth
- Short convalescence time
- Can be used on small, medium, large and giant breeds

### Key Learning Objectives

By the end of this course delegates should have an understanding of:




- Why we've developed another cruciate surgery technique
- Biomechanics and theoretical foundations of MMP
- Case selection for MMP
- How to perform the MMP procedure








\*As of October 2017



**Order by Phone**

 +44 (0) 845 045 0259  
 +1 772-562-6044  
 +61 (0) 8 9590 8850

 +49 (0) 4191 8030013  
 +34 902 47 47 01  
 +39 (0) 391 11 60 784

 +27 (0) 83 227 8181  
 +886 4 2258 6668

**Order Online**

All of our products are available to buy online at:

[www.orthomed.co.uk/nz](http://www.orthomed.co.uk/nz)

# RidgeStop™

## For Patella Luxation

### Course Content

RidgeStop™ is a novel surgical technique for treatment of patella luxation. A simple, affordable alternative to sulcoplasty surgery that can be used alone as the sole treatment for patella luxation, or as an adjunct to re-alignment operations.

The course is in two parts: a lecture based overview of the pathogenesis and treatment of patella luxation and a review of current surgical treatment options. This will be followed by an introduction to RidgeStop™, its application and case presentations, before a practical session where delegates will learn to use the RidgeStop™ system.

### Benefits

- Made from medical grade Ultra High Molecular Weight Polyethylene
- Minimal joint interference
- Medial or lateral application
- Minimally invasive
- Uses standard cortical screws





### Key Learning Objectives

By the end of this course delegates should have an understanding of:

- The aims of sulcoplasty
- Treatment selection in patella luxation
- The concept of RidgeStop™
- How to carry out the RidgeStop™ procedure





**Order by Phone**
 +44 (0) 845 045 0259

 +1 772-562-6044


 +61 (0) 8 9590 8850

 +49 (0) 4191 8030013

 +34 902 47 47 01

 +39 (0) 391 11 60 784

 +27 (0) 83 227 8181

 +886 4 2258 6668
**Order Online**

All of our products are available to buy online at:

[www.orthomed.co.uk/nz](http://www.orthomed.co.uk/nz)

# SOP™

## For Internal Fracture Repair

**Course Content**

The SOP™ locking plate system was developed by veterinary surgeons for use in veterinary orthopaedics. Locking plates are powerful and versatile implants but the mechanics and methods of application are significantly different from conventional plate systems. The course will cover locking plate mechanics, case selection, methods of application as well as modes of failure and a review of the extensive published literature on the SOP™ system.

**Benefits**

- Over 2,000 users worldwide\*
- Over 100 published papers
- More than 50,000 cases to-date\*
- Locks using standard cortical screws
- Does not require perfect contouring

**Key Learning Objectives**




By the end of this course delegates should have an understanding of:




- Features of locking plates.
- Case selection and clinical application of SOP™ system
- How to use SOP™ in a range of orthopaedic applications





\*As of October 2017

**Order by Phone**

 +44 (0) 845 045 0259  
 +1 772-562-6044  
 +61 (0) 8 9590 8850

 +49 (0) 4191 8030013  
 +34 902 47 47 01  
 +39 (0) 391 11 60 784

 +27 (0) 83 227 8181  
 +886 4 2258 6668

**Order Online**

All of our products are available to buy online at:

[www.orthomed.co.uk/nz](http://www.orthomed.co.uk/nz)

# Advanced Spinal Course

## Course Content

This three day conference will focus on advanced neurosurgical procedures of the cervical and thoracolumbar spine. It will consist of lectures including instruments, anaesthetic protocols, triage and step-by-step tricks and tips to six of the most common spinal procedures – T-L laminectomy, cervical slot decompression, spinal fracture stabilisation, atlantoaxial stabilisation, L7 – S1 decompression and stabilisation.


Attendees will choose and perform three laboratories from five offered. Cadavers will be available for all laboratories. The SOP™ locking plate system will be used for stabilisations. The attendees should already have basic knowledge of these surgical conditions, and basic surgical proficiency.

## Key Learning Objectives

By the end of this course delegates should have an understanding of:

- Advanced neurosurgical procedures of the cervical and thoracolumbar spine
- Use of instruments anaesthetic and analgesic protocols
- Step-by-step tricks and tips of the six most common spinal procedures

Two-day wet laboratory on cadavers to encompass these surgeries





Order by Phone

+44 (0) 845 045 0259  
+1 772-562-6044  
+61 (0) 8 9590 8850

+49 (0) 4191 8030013  
+34 902 47 47 01  
+39 (0) 391 11 60 784

+27 (0) 83 227 8181  
+886 4 2258 6668

Order Online

All of our products are available to buy online at:

[www.orthomed.co.uk/nz](http://www.orthomed.co.uk/nz)

*" Just wanted to say a huge thank you for today. Not only was it a great MMP and RidgeStop™ course but a really interesting discussion on the subject of cruciate “repair” and meniscal injury. What an amazing day."*

Heidi Birkinshaw



“ Although there is an increased choice within the veterinary orthopaedic industry, I choose Orthomed for many of my instruments, power tools and specialist equipment. The quality of engineering and the life expectancy of their tools gives me great confidence in surgery when using them. ”

**Robin Stratton**  
Accrington Vets4Pets





“ *Having worked with Orthomed for the last three years, I have found them a pleasure to work with. They're supportive, friendly, have a 'happy to help' attitude and their product range is of the highest quality.* ”

**Ben Orton**

Black Notley Veterinary Surgery







# MMP

## Canine Cruciate Repair System


01




User Guide available as  
download or hard copy

MMP is a simple procedure which requires minimal instruments and removes the need for void-filling or bone-grafting. MMP uses precision instruments combined with OrthoFoam™, a unique titanium foam wedge with an appropriately rigid, porous structure which allows new bone cells to proliferate throughout the wedge, bonding it to the host bone and creating a solid repair which considerably reduces convalescence and considerably improves the long term prognosis for the patient.

### In This Section






MMP Implants








MMP Instruments

20

**Order by Phone**

 +44 (0) 845 045 0259  
 +1 772-562-6044  
 +61 (0) 8 9590 8850

 +49 (0) 4191 8030013  
 +34 902 47 47 01  
 +39 (0) 391 11 60 784

 +27 (0) 83 227 8181  
 +886 4 2258 6668

**Order Online**

All of our products are available to buy online at:

[www.orthomed.co.uk/nz](http://www.orthomed.co.uk/nz)

# MMP

## Canine Cruciate Repair System

### About MMP

The Modified Maquet Procedure (MMP) is so named to acknowledge the instigator of these operations, Dr P Maquet, a human orthopaedic surgeon working in Belgium during the 1960s. MMP uses a wedge shaped implant of Titanium OrthoFoam™ which both defines the degree of advancement of the tibial tuberosity and holds the bone in its new place while the bony ingrowth that provides permanent, biomechanically robust fixation develops.

The Titanium OrthoFoam™ wedge provides a robust early fixation without the need for support bandages or a lengthy period of rest. The potentially disruptive forces that act to displace the distal end of the tibial tuberosity cranially following advancement are controlled using either a wire tension band, a titanium staple or a plate and screws removing the need for bone grafting or similar.



### About Titanium OrthoFoam™


The titanium foam is biocompatible and possesses mechanical properties which make it suitable for implantation in loadbearing situations.


The interconnected pores permit rapid bone ingrowth and fixation to surrounding bone. Orthomed began work with the National Research Council of Canada in 2008 which has resulted in the first application of OrthoFoam™ in the MMP wedge for canine cruciate disease.

The open porous structure and sympathetic mechanical characteristics of the Titanium OrthoFoam™ implant encourages early and sustained ingrowth of bone removing the need for bone grafting or similar.




**Order by Phone**

 +44 (0) 845 045 0259

 +1 772-562-6044

 +61 (0) 8 9590 8850

 +49 (0) 4191 8030013

 +34 902 47 47 01

 +39 (0) 391 11 60 784

 +27 (0) 83 227 8181

 +886 4 2258 6668

**Order Online**

All of our products are available to buy online at:

[www.orthomed.co.uk/nz](http://www.orthomed.co.uk/nz)

21

*"Since adopting the MMP as my preferred mode of treatment for ACL deficiency I have found all the products to be innovative, practical and very user friendly. Excellent customer service."*

Grahame Tasker  
Parkhill Vets

**Benefits of MMP****Simplify the surgical technique**

- Surgical approach
- Instruments
- Fixation
- Eliminates bone grafting/void-filling

**Shorten the convalescence**




- Robust fixation system
- "Less invasive" surgical technique




**Universal application**



- Appropriate for dogs of any size

**Also Available** - Animation, video, user guide, after-care guide

**Order by Phone**

 +44 (0) 845 045 0259  
 +1 772-562-6044  
 +61 (0) 8 9590 8850

 +49 (0) 4191 8030013  
 +34 902 47 47 01  
 +39 (0) 391 11 60 784

 +27 (0) 83 227 8181  
 +886 4 2258 6668

**Order Online**

All of our products are available to buy online at:

[www.orthomed.co.uk/nz](http://www.orthomed.co.uk/nz)

# MMP Resources



## An Owners Guide to MMP

We have produced an information booklet that allows the patient owner to understand the MMP operation; the benefits, what the procedure involves but also the important part that the owners play in giving the dog the best chance of an uneventful return to normal activity. These booklets are available by contacting your office.



## Staple Inserter

Available for both standard and small breed size staples we have developed a tool to aid staple insertion. This disposable tool is simple to use, holding the staple so that it can be tapped in until fully engaged.

## Papers



The following clinical MMP papers are available online and in print








***A review of 131 dogs with lameness due to CCL failure treated with MMP***

*Ness, Midgely and Harms (2012) Veterinary Orthopaedic Society Annual Conference, Crested Butte, CO, USA*

**Order by Phone**

 +44 (0) 845 045 0259  
 +1 772-562-6044  
 +61 (0) 8 9590 8850

 +49 (0) 4191 8030013  
 +34 902 47 47 01  
 +39 (0) 391 11 60 784

 +27 (0) 83 227 8181  
 +886 4 2258 6668

**Order Online**

All of our products are available to buy online at:


[www.orthomed.co.uk/nz](http://www.orthomed.co.uk/nz)

23

*“Orthomed is a great company with amazing products. Their product support system is top notch and they answer any questions you have quickly. I would definitely recommend Orthomed to anyone that has an interest in orthopaedics.”*

Carl Myers

Theodore Veterinary Hospital

**Staple Placement Guide**

An information sheet to visualise the positions for staple placement.

**Surgical Animations**

Visual walk-through of our products and how they are used during surgery.



**Planned wedge size compared to achieved advancement in dogs undergoing the modified Maquet procedure**

Kapler, Marcellin-Little and Roe (2015). VCOT 2015;28(6):379-384



**The Modified Maquet Procedure (MMP) in Dogs: Technical Development and Initial Clinical Experience**


Ness (2016). JAAHA 2016;52(4):242-250

# MMP Implants

Titanium OrthoFoam™ is suitable for dogs of any size, including large dogs, where historically the success of cruciate surgery has been limited. To date, over 35,000 procedures have been carried out and their success underpins the growing popularity of this novel procedure.









Resources available online.







## Order by Phone

 +44 (0) 845 045 0259  
 +1 772-562-6044  
 +61 (0) 8 9590 8850

 +49 (0) 4191 8030013  
 +34 902 47 47 01  
 +39 (0) 391 11 60 784

 +27 (0) 83 227 8181  
 +86 4 2258 6668

## Order Online

All of our products are available to buy online at:

[www.orthomed.co.uk/nz](http://www.orthomed.co.uk/nz)

25

# OrthoFoam™ MMP Implants



## Standard

Size	Code	Price
6mm	MMP054006	\$285.00
7.5mm	MMP0540075	\$285.00
9mm	MMP054009	\$285.00
10.5mm	MMP0540105	\$285.00

## Small Breed

Size	Code	Price
12mm	MMP054012	\$285.00
13.5mm	MMP0540135	\$285.00
15mm	MMP054015	\$285.00

Size	Code	Price
3mm	MMP154003	\$285.00
4mm	MMP154004	\$285.00
5mm	MMP154005	\$285.00
6mm	MMP154006	\$285.00

Size	Code	Price
7mm	MMP0154007	\$285.00
8mm	MMP0154008	\$285.00

# MMP Ti Staples and Inserters

## Standard

Size	Code	Price
6mm Leg	STAPLE0620TI	\$86.00
7mm Leg	STAPLE0720TI	\$86.00
Standard Inserter	STAPLEINS20	\$4.50



## Small Breed

Size	Code	Price
5mm Leg	STAPLE0513TI	\$86.00
6mm Leg	STAPLE0613TI	\$86.00
Small Breed Inserter	STAPLEINS13	\$4.50




MMP

# Instruments




The use of a carefully engineered saw guide gives control and ensures an osteotomy of the correct length, precisely directed and positioned to ensure an appropriately “thick” tibial tuberosity while at the same time protecting adjacent structures from iatrogenic injury.








Resources available online.



## Order by Phone

 +44 (0) 845 045 0259  
 +1 772-562-6044  
 +61 (0) 8 9590 8850

 +49 (0) 4191 8030013  
 +34 902 47 47 01  
 +39 (0) 391 11 60 784

 +27 (0) 83 227 8181  
 +886 4 2258 6668

## Order Online

All of our products are available to buy online at:

[www.orthomed.co.uk/nz](http://www.orthomed.co.uk/nz)

27

# MMP Instruments



## Standard

Product	Code	Price
<b>MMP Instrument Kit</b>	<b>MMP010000KIT</b>	<b>\$4,470.00</b>

MMP Instrument Kit contains the following

Product	Code	Price
<b>MMP Instrument Tray</b>	<b>MMPTRAY-STD</b>	<b>\$800.00</b>
<b>MMP Saw Guide</b> (incl. 1 proximal pin & 3 distal pins)	<b>MMP024100</b>	<b>\$490.00</b>
<b>MMP Wedge Inserter</b> (incl. long and short drill guide)	<b>MMP034100</b>	<b>\$1,190.00</b>
<b>1.5mm Drill Bit (JC)</b>	<b>900000015</b>	<b>\$40.00</b>
<b>2.0mm Drill Bit (JC)</b>	<b>900000020</b>	<b>\$40.00</b>
<b>3.5mm Short Drill Bit for MMP Saw Guide</b>	<b>MMP033001</b>	<b>\$60.00</b>
<b>MMP Combined K-Wire Bender and Staple Drill Guide</b>	<b>MMP010002-SDG</b>	<b>\$290.00</b>
<b>Staple System 2mm Pin</b>	<b>MMP010020</b>	<b>\$60.00</b>
<b>6mm MMP Trial Wedge</b>	<b>MMP044006</b>	<b>\$300.00</b>
<b>7.5mm MMP Trial Wedge</b>	<b>MMP044075</b>	<b>\$300.00</b>
<b>9mm MMP Trial Wedge</b>	<b>MMP044009</b>	<b>\$300.00</b>
<b>10.5mm MMP Trial Wedge</b>	<b>MMP044105</b>	<b>\$300.00</b>
<b>12mm MMP Trial Wedge</b>	<b>MMP044012</b>	<b>\$300.00</b>

Available as an add on:

Size	Code	Price
<b>13.5mm MMP Trial Wedge</b>	<b>MMP044135</b>	<b>\$300.00</b>
<b>15mm MMP Trial Wedge</b>	<b>MMP044015</b>	<b>\$300.00</b>

## Small Breed

Size	Code	Price
<b>MMP Small Breed Instrument Kit</b>	<b>MMP090000KIT</b>	<b>\$4,010.00</b>

MMP Small Breed Instrument Kit contains the following

Product	Code	Price
<b>MMP SB Instrument Tray</b>	<b>MMPTRAY-SB</b>	<b>\$700.00</b>
<b>MMP-SB Saw Guide</b> (incl. positioning pins & drill sleeve)	<b>MMP092310</b>	<b>\$600.00</b>
<b>MMP-SB Wedge Inserter</b> (incl. drill guide)	<b>MMP093110</b>	<b>\$990.00</b>
<b>1.5mm Drill Bit (JC)</b>	<b>900000015</b>	<b>\$40.00</b>
<b>2.0mm Drill Bit (JC)</b>	<b>900000020</b>	<b>\$40.00</b>
<b>2.5mm Short Drill Bit for MMP-SB Saw Guide</b>	<b>MMP133001</b>	<b>\$60.00</b>
<b>3.5mm Short Drill Bit for MMP Saw Guide</b>	<b>MMP033001</b>	<b>\$60.00</b>
<b>MMP Combined K-Wire Bender and Staple Drill Guide</b>	<b>MMP010002-SDG</b>	<b>\$290.00</b>
<b>Staple System 2mm Pin</b>	<b>MMP010020</b>	<b>\$60.00</b>
<b>3/4mm MMP-SB Trial Wedge</b>	<b>MMP094034</b>	<b>\$390.00</b>
<b>5/6mm MMP-SB Trial Wedge</b>	<b>MMP094056</b>	<b>\$390.00</b>
<b>7/8mm MMP-SB Trial Wedge</b>	<b>MMP094078</b>	<b>\$390.00</b>







# RidgeStop™

## For Patella Luxation


02




User Guide available as  
download or hard copy

RidgeStop™ is a radical new product which provides surgeons with a cost effective alternative to sulcoplasty when carrying out procedures used to correct patella luxation. It uses a synthetic rim, held in by standard cortical screws, when carrying out procedures to manage patella luxation. The synthetic ridge is secured with cortical bone screws and avoids further damage to the articular cartilage of the trochlear groove. It can be used as a stand-alone procedure in less severe cases, or applied following realignment surgeries such as distal femoral osteotomy or tibial crest translocation.

### In This Section



RidgeStop™ Implants



RidgeStop™ Instruments

30

**Order by Phone**

🇬🇧 +44 (0) 845 045 0259  
 🇺🇸 +1 772-562-6044  
 🇦🇺 +61 (0) 8 9590 8850

🇩🇪 +49 (0) 4191 8030013  
 🇪🇸 +34 902 47 47 01  
 🇮🇹 +39 (0) 391 11 60 784

🇳🇦 +27 (0) 83 227 8181  
 🇵🇰 +886 4 2258 6668

**Order Online**

All of our products are available to buy online at:

[www.orthomed.co.uk/nz](http://www.orthomed.co.uk/nz)

# RidgeStop™

## For Patella Luxation

### About RidgeStop™

RidgeStop™ can be used as a “stand-alone” procedure for less severe cases of patella luxation – Grade 2 or better. Alternatively, RidgeStop™ is employed as an adjunct procedure when re-alignment operations such as distal femoral osteotomy or tibial crest translocation have been performed. In such cases, RidgeStop™ acts effectively as an alternative to sulcoplasty procedures. When used as an adjunct to a re-alignment operation, the re-alignment surgery is completed before the RidgeStop™ is applied as part of the same anaesthetic procedure.

### Features

#### Long proximal arm

– Effective with patella alta

#### Deep axial geometry

– Optimal patella stability

#### Thinner abaxial edge

– Minimises joint interference

#### Asymmetry

– Medial or lateral ridge application

**UHMWPE (Ultra-high-molecular-weight polyethylene)**

**Used extensively in veterinary and human joint replacements**

**Standard cortical bone screws**

**Complex 3 dimensional geometry**

### Benefits of RidgeStop™

- ✓ Simple and quick
- ✓ Minimally traumatic
- ✓ “Stand-alone” procedure often possible, simplifying surgery
- ✓ Effective practical add on in cruciate surgeries where MPL needs to be addressed
- ✓ Reversible procedure
- ✓ Doesn't absolve surgeon from assessing need for re-alignment osteotomies

**Also Available** - Animation, video, user guide, after-care guide

**Order by Phone**

🇬🇧 +44 (0) 845 045 0259  
 🇺🇸 +1 772-562-6044  
 🇦🇺 +61 (0) 8 9590 8850

🇩🇪 +49 (0) 4191 8030013  
 🇪🇸 +34 902 47 47 01  
 🇮🇹 +39 (0) 391 11 60 784

🇿🇦 +27 (0) 83 227 8181  
 🇹🇼 +886 4 2258 6668

**Order Online**

All of our products are available to buy online at:

[www.orthomed.co.uk/nz](http://www.orthomed.co.uk/nz)

31

**Practical Tips**

- 💡 Medial arthrotomy approach
- 💡 No need to close joint capsule
- 💡 Make sure it does not cross patella line of travel
- 💡 Evaluate tracking before and after application
- 💡 In case of “stand-alone” application, always warn owners that further realignment procedure may still be needed



# RidgeStop™ Implants

The RidgeStop™ implants are precision manufactured using Ultra High Molecular-Weight Polyethylene (UHMWP) – the same material used extensively in veterinary and human joint replacement. This non-toxic, ultra-tough material transfers loads effectively and has an extremely low coefficient of friction making it ideal for patella luxation surgery.





Resources available online.






## Order by Phone

 +44 (0) 845 045 0259

 +1 772-562-6044


 +61 (0) 8 9590 8850

 +49 (0) 4191 8030013

 +34 902 47 47 01

 +39 (0) 391 11 60 784

 +27 (0) 83 227 8181

 +886 4 2258 6668

## Order Online

All of our products are available to buy online at:

[www.orthomed.co.uk/nz](http://www.orthomed.co.uk/nz)

33

# RidgeStop™ Implants



## LL RM Left Lateral, Right Medial

Size	Code	Price
1.5	RS15LLRM-LP (low profile)	\$321.00
1.5	RS15LLRM	\$321.00
2.0	RS20LLRM	\$330.00
2.5	RS25LLRM	\$340.00
3.0	RS30LLRM	\$349.00
3.5	RS35LLRM	\$359.00

## LM RL Left Medial, Right Lateral

Size	Code	Price
1.5	RS15LMRL-LP (low profile)	\$321.00
1.5	RS15LMRL	\$321.00
2.0	RS20LMRL	\$330.00
2.5	RS25LMRL	\$340.00
3.0	RS30LMRL	\$349.00
3.5	RS35LMRL	\$359.00


RidgeStop™

# Instruments

Each implant size has a corresponding drill guide to ensure ideal placement of the implant. Dependent on the implant size, they will require either two or three standard cortical bone screws. This guide allows not only ideal placement of the implant but also alignment of the divergent screw holes.



Resources available online.



## Order by Phone

UK +44 (0) 845 045 0259  
 USA +1 772-562-6044  
 AUS +61 (0) 8 9590 8850

GER +49 (0) 4191 8030013  
 ESP +34 902 47 47 01  
 ITA +39 (0) 391 11 60 784

SZA +27 (0) 83 227 8181  
 CHN +86 4 2258 6668

## Order Online

All of our products are available to buy online at:

[www.orthomed.co.uk/nz](http://www.orthomed.co.uk/nz)

35

# RidgeStop™ Instruments



## Starter Kit

Item	Code	Price
Starter kit for sizes 1.5, 2.0 and 2.5. Includes tray and six corresponding implants	RSKIT	\$5,312.00

## RidgeStop™ Tray

Size	Code	Price
RidgeStop™ Tray 1	RSTRAY1	\$614.00

## RidgeStop™ Pin

Size	Code	Price
RidgeStop™ Pin	RSDGPIN	\$28.00

## RidgeStop™ Drill Guides

**LL RM** Left Lateral, Right Medial

Size	Code	Price
1.5	RSDG15LLRM	\$415.00
2.0	RSDG20LLRM	\$425.00
2.5	RSDG25LLRM	\$434.00
3.0	RSDG30LLRM	\$444.00
3.5	RSDG35LLRM	\$453.00

**LM RL** Left Medial, Right Lateral

Size	Code	Price
1.5	RSDG15LMRL	\$415.00
2.0	RSDG20LMRL	\$425.00
2.5	RSDG25LMRL	\$434.00
3.0	RSDG30LMRL	\$444.00
3.5	RSDG35LMRL	\$453.00









# SOP™

## For Internal Fracture Repair


03




User Guide available as  
download or hard copy

DCPs and reconstruction plates are familiar to most orthopaedic surgeons, but in recent years, the emergence of locking plates has allowed for greater flexibility and more successful surgical outcomes. We have taken the technique one step further with the SOP (String of Pearls) System, which provides exceptionally flexible options to the surgeon, particularly in complicated situations where contouring on multiple planes is required or where the location of the fracture makes exact contouring difficult. The SOP™ system is the only locking plate system on the market that locks with conventional cortical bone screws.


### In This Section



SOP™ Locking Plates









SOP™ Instruments





SOP™ TPLO & TPO

**Order by Phone**

 +44 (0) 845 045 0259  
 +1 772-562-6044  
 +61 (0) 8 9590 8850

 +49 (0) 4191 8030013  
 +34 902 47 47 01  
 +39 (0) 391 11 60 784

 +27 (0) 83 227 8181  
 +886 4 2258 6668

**Order Online**

All of our products are available to buy online at:

[www.orthomed.co.uk/nz](http://www.orthomed.co.uk/nz)

# SOP™

The world leading locking plate system that uses standard cortical screws

**About SOP™**

The SOP™ (String of Pearls) system is a strong versatile locking plate system. The SOP™ plate consists of a series of cylindrical sections ('internodes') and spherical components ('pearls') that can be bent and twisted to suit almost every situation. SOP™ provides exceptionally flexible options to the surgeon, particularly in complicated situations where contouring on multiple planes is required or where the location of the fracture makes exact contouring difficult. As well as fracture repair applications also include ostetomies, spinal fusion, arthrodesis.

**Features****Available in 3 sizes:**

– 2.0mm, 2.7mm, 3.5mm

**Locks using conventional cortical bone screws**

**No need for lots of expensive equipment**

**Plates don't need to be contoured to exact position**

**Works as an internal external plate**



- ✓ 2,000 users worldwide\*
- ✓ Approximately 100 published papers in veterinary peer reviewed journals\*

- ✓ Over 100,000 cases to date\*
- ✓ Used by general practitioners, universities and referral centres

\*As of October 2017

**Order by Phone**

🇬🇧 +44 (0) 845 045 0259

🇺🇸 +1 772-562-6044

🇦🇺 +61 (0) 8 9590 8850

🇩🇪 +49 (0) 4191 8030013

🇪🇸 +34 902 47 47 01

🇮🇹 +39 (0) 391 11 60 784

🇿🇦 +27 (0) 83 227 8181

🇵🇹 +353 4 2258 6668

**Order Online**


All of our products are available to buy online at:

[www.orthomed.co.uk/nz](http://www.orthomed.co.uk/nz)

39

*"I have been using Orthomed SOP™ plates for eight years. They are an excellent adjunct to dynamic plate systems and allow for locking applications with standard cortical screws."*

Dr. Dan Degner

**Range of Applications**

- ✔ Fractures of the Femur
- ✔ Fractures of the Humerus
- ✔ Fractures of the Tibia
- ✔ Fractures of the Antebrachium
- ✔ Pelvic Fractures
- ✔ Plate - Rod Constructs
- ✔ Spinal Fractures
- ✔ Surgical Spine Distraction-Fusion

**Also Available** - Animation, video, user guide, after-care guide



SOP™

# Locking Plates


SOP™ Locking Plates	41
SOP™ Titanium Locking Plates	42

The SOP™ locking plate system can be contoured with 6 degrees of freedom and uses conventional screws. The plates come in three sizes – 2.0mm, 2.7mm and 3.5mm. Used in conjunction with the SOP™ instrumentation, the SOP™ locking plates can be easily contoured for a multitude of applications. With over 100,000 cases, and over 100 peer reviewed publications the SOP™ system is still the original and best.

The SOP™ titanium locking plate system works in the same way as the steel version, but uses titanium screws. Used for when MRI imaging is required.



Resources available online.







## Order by Phone


 +44 (0) 845 045 0259

 +49 (0) 4191 8030013


 +27 (0) 83 227 8181

 +1 772-562-6044

 +34 902 47 47 01

 +886 4 2258 6668

 +61 (0) 8 9590 8850

 +39 (0) 391 11 60 784

## Order Online

All of our products are available to buy online at:

[www.orthomed.co.uk/nz](http://www.orthomed.co.uk/nz)


41

# SOP™ Locking Plates

2.0mm



2.0mm Holes

6° Of Freedom  
Use Bending Irons (See Page 46)

Size	Code	Price
3 hole	800020003	\$50.00
4 hole	800020004	\$67.00
5 hole	800020005	\$82.00


Size	Code	Price
6 hole	800020006	\$99.00
7 hole	800020007	\$116.00
8 hole	800020008	\$133.00

Size	Code	Price
10 hole	800020010	\$165.00
12 hole	800020012	\$199.00
20 hole	800020020	\$331.00

2.7mm



2.7mm Holes

6° Of Freedom  
Use Bending Irons (See Page 46)Use Bending Tees  
(See Page 46)

Size	Code	Price
3 hole	800027003	\$50.00
4 hole	800027004	\$67.00
5 hole	800027005	\$82.00
6 hole	800027006	\$99.00

Size	Code	Price
7 hole	800027007	\$116.00
8 hole	800027008	\$133.00
9 hole	800027009	\$149.00
10 hole	800027010	\$165.00


Size	Code	Price
11 hole	800027011	\$182.00
12 hole	800027012	\$199.00
13 hole	800027013	\$215.00
14 hole	800027014	\$232.00

Size	Code	Price
15 hole	800027015	\$248.00
16 hole	800027016	\$265.00
20 hole	800027020	\$331.00

3.5mm



3.5mm Holes

6° Of Freedom  
Use Bending Irons (See Page 46)Use Bending Tees  
(See Page 46)




Size	Code	Price
3 hole	800035003	\$50.00
4 hole	800035004	\$67.00
5 hole	800035005	\$82.00
6 hole	800035006	\$99.00
7 hole	800035007	\$116.00




Size	Code	Price
8 hole	800035008	\$133.00
9 hole	800035009	\$149.00
10 hole	800035010	\$165.00
11 hole	800035011	\$182.00
12 hole	800035012	\$199.00



Size	Code	Price
13 hole	800035013	\$215.00
14 hole	800035014	\$232.00
15 hole	800035015	\$248.00
16 hole	800035016	\$265.00
18 hole	800035018	\$298.00

Size	Code	Price
20 hole	800035020	\$331.00
30 hole	800035030	\$496.00

**Order by Phone**

 +44 (0) 845 045 0259  
 +1 772-562-6044  
 +61 (0) 8 9590 8850

 +49 (0) 4191 8030013  
 +34 902 47 47 01  
 +39 (0) 391 11 60 784

 +27 (0) 83 227 8181  
 +886 4 2258 6668

**Order Online**

All of our products are available to buy online at:


[www.orthomed.co.uk/nz](http://www.orthomed.co.uk/nz)

## SOP™ Titanium Locking Plates

**2.0mm**



2.7mm Holes



6° Of Freedom  
Use Bending Irons (See Page 46)



Size	Code	Price
3 hole	T800020003	\$61.00
4 hole	T800020004	\$81.00

Size	Code	Price
5 hole	T800020005	\$102.00
6 hole	T800020006	\$122.00

Size	Code	Price
7 hole	T800020007	\$143.00
8 hole	T800020008	\$163.00

Size	Code	Price
10 hole	T800020010	\$204.00
12 hole	T800020012	\$245.00

**2.7mm**



2.7mm Holes



6° Of Freedom  
Use Bending Irons (See Page 46)



Use Ti Bending Tees  
(See Page 46)



Size	Code	Price
3 hole	T800027003	\$61.00
4 hole	T800027004	\$81.00
5 hole	T800027005	\$102.00
6 hole	T800027006	\$122.00

Size	Code	Price
7 hole	T800027007	\$143.00
8 hole	T800027008	\$163.00
9 hole	T800027009	\$183.00
10 hole	T800027010	\$204.00

Size	Code	Price
11 hole	T800027011	\$224.00
12 hole	T800027012	\$245.00
13 hole	T800027013	\$265.00
14 hole	T800027014	\$286.00

Size	Code	Price
15 hole	T800027015	\$306.00
16 hole	T800027016	\$327.00

**3.5mm**



3.5mm Holes



6° Of Freedom  
Use Bending Irons (See Page 46)



Use Ti Bending Tees  
(See Page 46)



Size	Code	Price
3 hole	T800035003	\$61.00
4 hole	T800035004	\$81.00
5 hole	T800035005	\$102.00
6 hole	T800035006	\$122.00

Size	Code	Price
7 hole	T800035007	\$143.00
8 hole	T800035008	\$163.00
9 hole	T800035009	\$183.00
10 hole	T800035010	\$204.00

Size	Code	Price
11 hole	T800035011	\$224.00
12 hole	T800035012	\$245.00
13 hole	T800035013	\$265.00
14 hole	T800035014	\$286.00

Size	Code	Price
15 hole	T800035015	\$306.00
16 hole	T800035016	\$327.00
18 hole	T800035018	\$367.00
20 hole	T800035020	\$408.00



SOP™


# Instruments

SOP™ Starter Kits	45
SOP™ Bending Irons	46
SOP™ Drill Guides	46
SOP™ Bending Tees	46
SOP™ Bending Tees Case	47
SOP™ Plate Holder	47
SOP™ Instrument Tray	47
SOP™ Wire Guides	47



Resources available online.

Compromising of drill guide, bending irons and tees, this is all the instrumentation you need to get started with using the SOP™ system. They can be bought separately, or as a starter kit for each size.





## Order by Phone

🇬🇧 +44 (0) 845 045 0259

🇺🇸 +1 772-562-6044

🇦🇺 +61 (0) 8 9590 8850

🇩🇪 +49 (0) 4191 8030013

🇪🇸 +34 902 47 47 01

🇮🇹 +39 (0) 391 11 60 784

🇳🇬 +27 (0) 83 227 8181

🇿🇦 +886 4 2258 6668

## Order Online

All of our products are available to buy online at:

[www.orthomed.co.uk/nz](http://www.orthomed.co.uk/nz)

45

# SOP™ Starter Kits



2.0mm



Size	Code	Price
2.0mm	800000020	\$1,239.00
2.0mm (with tray)	800000020TRAY	\$1,889.00

12 hole plate included

2.7mm



Size	Code	Price
2.7mm	800000027	\$1,463.00
2.7mm (with tray)	800000027TRAY	\$2,213.00

16 hole plate, bending tees, & tees case included




3.5mm








Size	Code	Price
3.5mm	800000035	\$1,529.00
3.5mm (with tray)	800000035TRAY	\$2,354.00

20 hole plate, bending tees, & tees case included

## Order by Phone

 +44 (0) 845 045 0259  
 +1 772-562-6044  
 +61 (0) 8 9590 8850

 +49 (0) 4191 8030013  
 +34 902 47 47 01  
 +39 (0) 391 11 60 784

 +27 (0) 83 227 8181  
 +886 4 2258 6668

## Order Online

All of our products are available to buy online at:

[www.orthomed.co.uk/nz](http://www.orthomed.co.uk/nz)

## SOP™ Bending Irons



2.0mm

Size	Code	Price
Left	800400020	\$380.00
Right	800410020	\$380.00

2.7mm

Size	Code	Price
Left	800400027	\$380.00
Right	800410027	\$380.00

3.5mm

Size	Code	Price
Left	800400035	\$380.00
Right	800410035	\$380.00

## SOP™ Drill Guides

Size	Code	Price
1.5/2.0mm	800201520	\$280.00
2.0/2.7mm	800202027	\$280.00
2.5/3.5mm	800202535	\$280.00



## SOP™ Bending Tees



Use 2.5mm Hex Screwdriver  
(See Pages 77/78)






Standard




Size	Code	Price
2.7mm (pack of 20)	800400127	\$90.00
3.5mm (pack of 20)	800400135	\$90.00



Titanium

Size	Code	Price
2.7mm (pack of 20)	T800400127	\$93.00
3.5mm (pack of 20)	T800400135	\$93.00

## Order by Phone

 +44 (0) 845 045 0259  
 +1 772-562-6044  
 +61 (0) 8 9590 8850

 +49 (0) 4191 8030013  
 +34 902 47 47 01  
 +39 (0) 391 11 60 784

 +27 (0) 83 227 8181  
 +86 4 2258 6668

## Order Online

All of our products are available to buy online at:

[www.orthomed.co.uk/nz](http://www.orthomed.co.uk/nz)

47

## SOP™ Bending Tees Case

Size	Code	Price
2.7mm (includes insert tray)	800400101	\$68.00
3.5mm (includes insert tray)	800400100	\$68.00
Insert Tray	800400000	\$13.20



## SOP™ Plate Holder

Size	Code	Price
2.0, 2.7 and 3.5mm	800300000	\$575.00



## SOP™ Instrument Tray

Size	Code	Price
2.0mm	SOPTRAY-20	\$650.00
2.7mm	SOPTRAY-27	\$750.00
3.5mm	SOPTRAY-35	\$825.00



## SOP™ Wire Guides

Size	Code	Price
2.0mm (Uses 0.8mm k-wire)	800420080	\$28.00
2.7mm (Uses 1.0mm k-wire)	800427010	\$28.00
3.5mm (Uses 1.0mm k-wire)	800435010	\$28.00



SOP™


# TPLO & TPO

Our SOP™ TPLO and TPO plates use the same technology as the regular SOP™ plates.

The SOP™ TPLO locking plate is both anatomical and symmetrical, the proximal part is pre-contoured and requires no further contouring, effective screw orientation is designed into the implant, whilst the TPO plates also uses the proven technology of SOP™. Both require the use of the SOP™ drill guide.



Resources available online.





**Order by Phone**

UK +44 (0) 845 045 0259  
 USA +1 772-562-6044  
 AUS +61 (0) 8 9590 8850

GER +49 (0) 4191 8030013  
 ESP +34 902 47 47 01  
 IT +39 (0) 391 11 60 784

SOUTH AFRICA +27 (0) 83 227 8181  
 CHINA +86 4 2258 6668

**Order Online**

All of our products are available to buy online at:

[www.orthomed.co.uk/nz](http://www.orthomed.co.uk/nz)

49

# SOP™ TPLO



2.7/3.5mm Holes



Use Cortical Screws (See Page 58)

The laser marked dot should lie directly over the osteotomy to allow correct placement and alignment of the plate



2.7mm

3.5mm

3.5mm Broad

Size	Code	Price
2.7mm	801027001	\$134.00

Size	Code	Price
3.5mm	801035001	\$134.00

Size	Code	Price
3.5mm Broad	801035002	\$189.00

# SOP™ TPO 3.5mm



Use 3.5mm Screws (See Page 58)



20°

25°

30°

Side	Code	Price
Left	81020100	\$200.00
Right	81020200	\$200.00

Side	Code	Price
Left	81025100	\$200.00
Right	81025200	\$200.00

Side	Code	Price
Left	81030100	\$200.00
Right	81030200	\$200.00



# Carpal Arthrodesis

## For Canine Carpus

---


04



User Guide available as  
download or hard copy

These implants have been carefully designed to minimise complications. There is no need to bend the plate, which helps to maintain its integrity and the proximal portion is engineered to be shorter and wider to accommodate the screws, before thinning at the proximal end to minimise stress. Potential catastrophic complications arising from the use of a cast are removed as patients operated on using our implants do not require a cast

### In This Section




PCA & ParCA

---




# Carpal Arthrodesis **PCA & ParCA**

The castles partial arthrodesis plates are the first range of plates designed specifically for canine partial carpal arthrodesis. The plates have several custom design features which make them ideal for surgery.









Resources available online.







**Order by Phone**

 +44 (0) 845 045 0259  
 +1 772-562-6044  
 +61 (0) 8 9590 8850

 +49 (0) 4191 8030013  
 +34 902 47 47 01  
 +39 (0) 391 11 60 784

 +27 (0) 83 227 8181  
 +886 4 2258 6668

**Order Online**

All of our products are available to buy online at:

[www.orthomed.co.uk/nz](http://www.orthomed.co.uk/nz)

53

# Castless PCA

2.0/1.5mm



Use 2.0/1.5mm Screws  
(See Page 58)



Size	Code	Price
55mm	107511055	\$165.00
60mm	107511060	\$170.00

Size	Code	Price
65mm	107511065	\$175.00

2.7/2.0mm



Use 2.7/2.0mm Screws  
(See Page 58)



Size	Code	Price
60mm	107510060	\$145.00
70mm	107510070	\$145.00

Size	Code	Price
75mm	107510075	\$154.00
85mm	107510085	\$163.00

3.5/2.7mm



Use 3.5/2.7mm Screws  
(See Page 58)



Size	Code	Price
110mm	107500110	\$152.00
120mm	107500120	\$152.00

Size	Code	Price
130mm	107500130	\$161.00
140mm	107500140	\$171.00

# Castless ParCA



Use 2.0/2.4mm Screws  
(See Page 58)




Use 3.5/2.7mm Screws  
(See Page 58)

Size	Code	Price
38mm	107520038	\$171.00
43mm	107520043	\$171.00
48mm	107520048	\$171.00
60mm	107520060	\$207.00
70mm	107520070	\$207.00
80mm	107520080	\$226.00
90mm	107520090	\$226.00



**Suitable screws sizes:**

	107520038	107520043	107520048	107520060	107520070	107520080	107520090
2	2.0 or 2.4 or 2.7mm			2.4 or 2.7mm		2.7 or 3.5mm	
1 & 3	2.0mm			2.7mm		2.7mm	
4 - 9	2.0mm			2.7mm		2.7mm	



# TPLO


## Tibial Plate Levelling Osteotomy

---


# 05

Tibial Plate Levelling Osteotomy reduces pain and instability in the canine stifle by removing a curved section of bone from the proximal end of the tibia, and rotating the bone to reduce the angle of the slope. Our TPLO plates are available in 3 sizes; 2.7mm, 3.5mm and 3.5mm Broad all pre-contoured to the top of the tibia with minimal further contouring to achieve a perfect fit. The proximal hole configuration of the 3.5mm Broad plate has been updated to minimise the risk of crossing screws, which reduces surgery time and considerably simplifies the procedure. Our TPLO plates are precision engineered from medical-grade stainless steel and laser cut to ensure that the finished product meets and exceeds current standards.

### In This Section






TPLO Plates








SOP™ TPLO Plates

---

**Order by Phone**

 +44 (0) 845 045 0259  
 +1 772-562-6044  
 +61 (0) 8 9590 8850

 +49 (0) 4191 8030013  
 +34 902 47 47 01  
 +39 (0) 391 11 60 784

 +27 (0) 83 227 8181  
 +886 4 2258 6668

**Order Online**

All of our products are available to buy online at:

[www.orthomed.co.uk/nz](http://www.orthomed.co.uk/nz)

# TPLO Plates

## Standard

Size	Code	Price
2.7mm Left	107401012	\$92.00
2.7mm Right	107401013	\$92.00

Size	Code	Price
3.5mm Left	107401011	\$107.00
3.5mm Right	107401010	\$107.00

Size	Code	Price
3.5mm Broad Left	107401015	\$142.00
3.5mm Broad Right	107401014	\$142.00



## Pre-Bent

Size	Code	Price
2.7mm Pre-Bent Left	107401004	\$92.00
2.7mm Pre-Bent Right	107401005	\$92.00




Size	Code	Price
3.5mm Pre-Bent Left	107401008	\$107.00
3.5mm Pre-Bent Right	107401009	\$107.00




Size	Code	Price
3.5mm Pre-Bent Broad Left	107401006	\$142.00
3.5mm Pre-Bent Broad Right	107401007	\$142.00







**Order by Phone**

 +44 (0) 845 045 0259  
 +1 772-562-6044  
 +61 (0) 8 9590 8850

 +49 (0) 4191 8030013  
 +34 902 47 47 01  
 +39 (0) 391 11 60 784

 +27 (0) 83 227 8181  
 +886 4 2258 6668

**Order Online**

All of our products are available to buy online at:

[www.orthomed.co.uk/nz](http://www.orthomed.co.uk/nz)

57

# SOP™ TPLO



2.7/3.5mm  
Holes



Use Cortical Screws  
(See Page 58)




Userguide  
Available

Size	Code	Price
2.7mm	801027001	\$134.00
3.5mm	801035001	\$134.00
3.5mm Broad	801035002	\$189.00



The laser marked dot should lie directly over the osteotomy to allow correct placement and alignment of the plate





1.5  
2.0  
2.5  
2.7  
3.5

10  
12  
14  
16  
18  
20  
22  
24  
26  
28  
30  
32  
34  
36  
38  
40  
45  
50  
55  
60



# Internal Fixation

---

## 06

We provide a comprehensive range of high quality screws, plates and washers. All products are manufactured from premium 316LVM (surgical grade) stainless steel or from titanium. The range includes broad and narrow compression plates, t-plates and acetabular plates, cortical nuts and washers, all manufactured to the highest quality, and in accordance with relevant ASTM and ISO standards.

### In This Section

Self Tapping Cortical Screws	60	Self Compression Plates	67
Standard Cortical Screws	62	T Plates	68
Titanium Self Tapping Cortical Screws	64	Acetabular Plates	69
Part Threaded Cancellous Screws	65	Cortical Nuts	69
Fully Threaded Cancellous Screws	65	Washers	69
Screw Racks	66	Cerclage Wire	70
Straight Plates	66	Steinmann Pins	70
Limited Contact Plates	66	Trochar Pins	71


---



60

## Order by Phone


 +44 (0) 845 045 0259

 +1 772-562-6044

 +61 (0) 8 9590 8850

 +49 (0) 4191 8030013

 +34 902 47 47 01

 +39 (0) 391 11 60 784

 +27 (0) 83 227 8181

 +886 4 2258 6668


## Order Online

All of our products are available to buy online at:

[www.orthomed.co.uk/nz](http://www.orthomed.co.uk/nz)

# Self Tapping Cortical Screws

1.5mm




Size	Code	Price
6mm	930515006	\$8.50
7mm	930515007	\$8.50
8mm	930515008	\$8.50
9mm	930515009	\$8.50

Size	Code	Price
10mm	930515010	\$9.00
11mm	930515011	\$9.00
12mm	930515012	\$9.00
13mm	930515013	\$9.00

Size	Code	Price
14mm	930515014	\$10.50
16mm	930515016	\$10.50
18mm	930515018	\$11.50
20mm	930515020	\$11.50

Size	Code	Price
22mm	930515022	\$11.50
24mm	930515024	\$12.00

2.0mm



Size	Code	Price
6mm	930520006	\$8.50
8mm	930520008	\$8.50
10mm	930520010	\$8.50
12mm	930520012	\$9.00
14mm	930520014	\$10.50

Size	Code	Price
16mm	930520016	\$10.50
18mm	930520018	\$10.50
20mm	930520020	\$11.50
22mm	930520022	\$11.50
24mm	930520024	\$12.50

Size	Code	Price
26mm	930520026	\$12.50
28mm	930520028	\$13.50
30mm	930520030	\$13.50
32mm	930520032	\$13.50
34mm	930520034	\$14.00

Size	Code	Price
36mm	930520036	\$14.00
38mm	930520038	\$15.00
40mm	930520040	\$15.00

2.4mm



Size	Code	Price
6mm	930524006	\$8.50
8mm	930524008	\$8.50
10mm	930524010	\$8.50
12mm	930524012	\$9.00
14mm	930524014	\$10.50

Size	Code	Price
16mm	930524016	\$10.50
18mm	930524018	\$10.50
20mm	930524020	\$11.50
22mm	930524022	\$11.50
24mm	930524024	\$12.50

Size	Code	Price
26mm	930524026	\$12.50
28mm	930524028	\$13.50
30mm	930524030	\$13.50
32mm	930524032	\$13.50
34mm	930524034	\$14.00

Size	Code	Price
36mm	930524036	\$14.00
38mm	930524038	\$15.00
40mm	930524040	\$15.00

2.7mm



Size	Code	Price
6mm	930527006	\$8.00
8mm	930527008	\$8.00

Size	Code	Price
10mm	930527010	\$8.00
12mm	930527012	\$8.00

Size	Code	Price
14mm	930527014	\$8.00
16mm	930527016	\$8.00

Size	Code	Price
18mm	930527018	\$9.00
20mm	930527020	\$9.00



**Order by Phone**

+44 (0) 845 045 0259  
 +1 772-562-6044  
 +61 (0) 8 9590 8850

+49 (0) 4191 8030013  
 +34 902 47 47 01  
 +39 (0) 391 11 60 784

+27 (0) 83 227 8181  
 +86 4 2258 6668

**Order Online**

All of our products are available to buy online at:

[www.orthomed.co.uk/nz](http://www.orthomed.co.uk/nz)

61

Size	Code	Price
22mm	930527022	\$9.00
24mm	930527024	\$9.00
26mm	930527026	\$9.00

Size	Code	Price
28mm	930527028	\$9.00
30mm	930527030	\$9.00
32mm	930527032	\$10.00

Size	Code	Price
34mm	930527034	\$10.00
36mm	930527036	\$10.00
38mm	930527038	\$10.00

Size	Code	Price
40mm	930527040	\$10.00

3.5mm



Size	Code	Price
10mm	930535010	\$9.00
12mm	930535012	\$9.00
14mm	930535014	\$9.00
16mm	930535016	\$9.00
18mm	930535018	\$10.00
20mm	930535020	\$10.00

Size	Code	Price
22mm	930535022	\$10.00
24mm	930535024	\$10.00
26mm	930535026	\$11.00
28mm	930535028	\$11.00
30mm	930535030	\$11.00
32mm	930535032	\$12.00

Size	Code	Price
34mm	930535034	\$12.00
36mm	930535036	\$12.00
38mm	930535038	\$12.00
40mm	930535040	\$12.00
42mm	930535042	\$13.00
44mm	930535044	\$13.00

Size	Code	Price
45mm	930535045	\$14.00
46mm	930535046	\$14.00
48mm	930535048	\$14.00
50mm	930535050	\$14.00
55mm	930535055	\$16.00
60mm	930535060	\$16.00

4.5mm



Size	Code	Price
12mm	930545012	\$10.50
14mm	930545014	\$10.50
16mm	930545016	\$10.50
18mm	930545018	\$11.50
20mm	930545020	\$11.50
22mm	930545022	\$11.50
24mm	930545024	\$11.50
26mm	930545026	\$11.50
28mm	930545028	\$12.00

Size	Code	Price
30mm	930545030	\$12.00
32mm	930545032	\$12.00
34mm	930545034	\$12.00
36mm	930545036	\$13.00
38mm	930545038	\$13.00
40mm	930545040	\$13.00
42mm	930545042	\$13.00
44mm	930545044	\$13.00
46mm	930545046	\$14.00

Size	Code	Price
48mm	930545048	\$14.00
50mm	930545050	\$14.00
52mm	930545052	\$14.00
54mm	930545054	\$14.00
55mm	930545055	\$15.00
56mm	930545056	\$15.00
58mm	930545058	\$16.00
60mm	930545060	\$17.00
65mm	930545065	\$18.00

Size	Code	Price
70mm	930545070	\$20.00
75mm	930545075	\$21.00
80mm	930545080	\$22.00
85mm	930545085	\$23.00
90mm	930545090	\$24.00
95mm	930545095	\$25.00
100mm	930545100	\$26.00
105mm	930545105	\$28.00
110mm	930545110	\$30.00

62

## Order by Phone

UK +44 (0) 845 045 0259  
 USA +1 772-562-6044  
 AUS +61 (0) 8 9590 8850

GER +49 (0) 4191 8030013  
 SPA +34 902 47 47 01  
 ITA +39 (0) 391 11 60 784

IND +27 (0) 83 227 8181  
 SAU +886 4 2258 6668


## Order Online

All of our products are available to buy online at:

[www.orthomed.co.uk/nz](http://www.orthomed.co.uk/nz)

# Standard Cortical Screws

1.5mm





Size	Code	Price
6mm	930215006	\$7.80
7mm	930215007	\$7.80
8mm	930215008	\$7.80
9mm	930215009	\$7.80

Size	Code	Price
10mm	930215010	\$8.20
11mm	930215011	\$8.20
12mm	930215012	\$8.20
13mm	930215013	\$8.20

Size	Code	Price
14mm	930215014	\$9.60
16mm	930215016	\$9.60
18mm	930215018	\$9.60
20mm	930215020	\$10.40

Size	Code	Price
22mm	930215022	\$10.40

2.0mm




Size	Code	Price
6mm	930220006	\$7.80
8mm	930220008	\$7.80
10mm	930220010	\$7.80
12mm	930220012	\$8.20

Size	Code	Price
14mm	930220014	\$9.60
16mm	930220016	\$9.60
18mm	930220018	\$9.60
20mm	930220020	\$10.40

Size	Code	Price
22mm	930220022	\$10.40
24mm	930220024	\$11.40
26mm	930220026	\$11.40
28mm	930220028	\$12.20

Size	Code	Price
30mm	930220030	\$12.20
32mm	930220032	\$12.20
34mm	930220034	\$12.60

2.7mm




Size	Code	Price
6mm	930227006	\$7.20
8mm	930227008	\$7.20
10mm	930227010	\$7.20
12mm	930227012	\$7.20
14mm	930227014	\$7.20

Size	Code	Price
16mm	930227016	\$7.20
18mm	930227018	\$8.20
20mm	930227020	\$8.20
22mm	930227022	\$8.20
24mm	930227024	\$8.20

Size	Code	Price
26mm	930227026	\$8.20
28mm	930227028	\$8.20
30mm	930227030	\$8.20
32mm	930227032	\$9.00
34mm	930227034	\$9.00

Size	Code	Price
36mm	930227036	\$9.00
38mm	930227038	\$9.00
40mm	930227040	\$9.00

3.5mm






Size	Code	Price
10mm	930235010	\$8.20




Size	Code	Price
12mm	930235012	\$8.20



Size	Code	Price
14mm	930235014	\$8.20

Size	Code	Price
16mm	930235016	\$8.20

**Order by Phone**

 +44 (0) 845 045 0259  
 +1 772-562-6044  
 +61 (0) 8 9590 8850

 +49 (0) 4191 8030013  
 +34 902 47 47 01  
 +39 (0) 391 11 60 784

 +27 (0) 83 227 8181  
 +886 4 2258 6668

**Order Online**

All of our products are available to buy online at:

[www.orthomed.co.uk/nz](http://www.orthomed.co.uk/nz)



Size	Code	Price
18mm	930235018	\$9.00
20mm	930235020	\$9.00
22mm	930235022	\$9.00
24mm	930235024	\$9.00
26mm	930235026	\$10.00

Size	Code	Price
28mm	930235028	\$10.00
30mm	930235030	\$10.00
32mm	930235032	\$10.80
34mm	930235034	\$10.80
36mm	930235036	\$10.80

Size	Code	Price
38mm	930235038	\$10.80
40mm	930235040	\$10.80
42mm	930235042	\$11.80
44mm	930235044	\$11.80
45mm	930235045	\$12.60

Size	Code	Price
46mm	930235046	\$12.60
48mm	930235048	\$12.60
50mm	930235050	\$12.60
55mm	930235055	\$14.40
60mm	930235060	\$14.40

**4.5mm**



Size	Code	Price
12mm	930245012	\$9.60
14mm	930245014	\$9.60
16mm	930245016	\$9.60
18mm	930245018	\$10.40
20mm	930245020	\$10.40
22mm	930245022	\$10.40
24mm	930245024	\$10.40
26mm	930245026	\$10.40
28mm	930245028	\$10.80

Size	Code	Price
30mm	930245030	\$10.80
32mm	930245032	\$10.80
34mm	930245034	\$10.80
36mm	930245036	\$11.80
38mm	930245038	\$11.80
40mm	930245040	\$11.80
42mm	930245042	\$11.80
44mm	930245044	\$11.80
46mm	930245046	\$12.60

Size	Code	Price
48mm	930245048	\$12.60
50mm	930245050	\$12.60
52mm	930245052	\$12.60
54mm	930245054	\$12.60
56mm	930245056	\$13.60
58mm	930245058	\$14.40
60mm	930245060	\$15.40
65mm	930245065	\$16.20
70mm	930245070	\$18.00

Size	Code	Price
75mm	930245075	\$19.00
80mm	930245080	\$19.80
85mm	930245085	\$20.80
90mm	930245090	\$21.60
95mm	930245095	\$22.60
100mm	930245100	\$23.40
105mm	930245105	\$25.20
110mm	930245110	\$27.00

**5.5mm**






Size	Code	Price
20mm	930255020	\$20.00
22mm	930255022	\$20.00
24mm	930255024	\$20.00
26mm	930255026	\$20.00
28mm	930255028	\$24.00
30mm	930255030	\$24.00
32mm	930255032	\$24.00
34mm	930255034	\$28.00
36mm	930255036	\$28.00




Size	Code	Price
38mm	930255038	\$28.00
40mm	930255040	\$32.00
42mm	930255042	\$32.00
44mm	930255044	\$32.00
46mm	930255046	\$36.00
48mm	930255048	\$36.00
50mm	930255050	\$40.00
52mm	930255052	\$40.00
54mm	930255054	\$40.00



Size	Code	Price
56mm	930255056	\$40.00
58mm	930255058	\$44.00
60mm	930255060	\$44.00
62mm	930255062	\$44.00
64mm	930255064	\$44.00
65mm	930255065	\$44.00
66mm	930255066	\$44.00
68mm	930255068	\$44.00
70mm	930255070	\$48.00

Size	Code	Price
72mm	930255072	\$48.00
74mm	930255074	\$48.00
76mm	930255075	\$48.00
78mm	930255076	\$48.00
80mm	930255080	\$52.00
85mm	930255085	\$52.00
90mm	930255090	\$56.00
95mm	930255095	\$56.00
10mm	930255100	\$60.00

**Order by Phone**

 +44 (0) 845 045 0259  
 +1 772-562-6044  
 +61 (0) 8 9590 8850

 +49 (0) 4191 8030013  
 +34 902 47 47 01  
 +39 (0) 391 11 60 784

 +27 (0) 83 227 8181  
 +886 4 2258 6668

**Order Online**

All of our products are available to buy online at:

[www.orthomed.co.uk/nz](http://www.orthomed.co.uk/nz)

## Titanium Self Tapping Cortical Screws

2.0mm



Size	Code	Price
6mm	T930520006	\$20.00
8mm	T930520008	\$20.00
10mm	T930520010	\$20.00

Size	Code	Price
12mm	T930520012	\$20.00
14mm	T930520014	\$20.00
16mm	T930520016	\$20.00

Size	Code	Price
18mm	T930520018	\$20.00
20mm	T930520020	\$20.00
22mm	T930520022	\$20.00

Size	Code	Price
24mm	T930520024	\$20.00

2.7mm



Size	Code	Price
8mm	T930527008	\$20.00
10mm	T930527010	\$20.00
12mm	T930527012	\$20.00
14mm	T930527014	\$20.00
16mm	T930527016	\$20.00

Size	Code	Price
18mm	T930527018	\$20.00
20mm	T930527020	\$20.00
22mm	T930527022	\$20.00
24mm	T930527024	\$20.00
26mm	T930527026	\$20.00

Size	Code	Price
28mm	T930527028	\$20.00
30mm	T930527030	\$20.00
32mm	T930527032	\$20.00
34mm	T930527034	\$20.00
36mm	T930527036	\$20.00

Size	Code	Price
38mm	T930527038	\$20.00
40mm	T930527040	\$20.00

3.5mm



Size	Code	Price
10mm	T930535010	\$20.00
12mm	T930535012	\$20.00
14mm	T930535014	\$20.00
16mm	T930535016	\$20.00
18mm	T930535018	\$20.00

Size	Code	Price
20mm	T930535020	\$20.00
22mm	T930535022	\$20.00
24mm	T930535024	\$20.00
26mm	T930535026	\$20.00
28mm	T930535028	\$20.00

Size	Code	Price
30mm	T930535030	\$20.00
32mm	T930535032	\$20.00
34mm	T930535034	\$20.00
36mm	T930535036	\$20.00
38mm	T930535038	\$20.00

Size	Code	Price
40mm	T930535040	\$20.00
45mm	T930535045	\$20.00
50mm	T930535050	\$20.00
55mm	T930535055	\$20.00
60mm	T930535060	\$20.00



**Order by Phone**

UK +44 (0) 845 045 0259  
 USA +1 772-562-6044  
 AUS +61 (0) 8 9590 8850

GER +49 (0) 4191 8030013  
 ESP +34 902 47 47 01  
 IT +39 (0) 391 11 60 784

SOUTH AFRICA +27 (0) 83 227 8181  
 CHINA +86 4 2258 6668


**Order Online**

All of our products are available to buy online at:


[www.orthomed.co.uk/nz](http://www.orthomed.co.uk/nz)

# Part Threaded Cancellous Screws


4.0mm




4.0mm  
Cancellous Screw




Standard  
Thread



2.5mm Hex Screwdriver  
(See Page 77/78)




2.5mm Pilot Hole  
(See Page 76)




Size	Code	Price	Size	Code	Price	Size	Code	Price	Size	Code	Price
10mm	931004010	\$8.00	20mm	931004020	\$9.00	30mm	931004030	\$11.00	55mm	931004055	\$14.00
12mm	931004012	\$8.00	22mm	931004022	\$9.00	35mm	931004035	\$11.00	60mm	931004060	\$14.00
14mm	931004014	\$8.50	24mm	931004024	\$10.00	40mm	931004040	\$11.00			
16mm	931004016	\$8.50	26mm	931004026	\$10.00	45mm	931004045	\$13.00			
18mm	931004018	\$9.00	28mm	931004028	\$10.00	50mm	931004050	\$13.00			

# Fully Threaded Cancellous Screws


4.0mm




4.0mm  
Fully Threaded Cancellous Screw




Standard  
Thread



2.5mm Hex Screwdriver  
(See Page 77/78)




2.5mm Pilot Hole  
(See Page 76)




Size	Code	Price	Size	Code	Price	Size	Code	Price	Size	Code	Price
10mm	931140010	\$7.00	20mm	931140020	\$8.50	30mm	931140030	\$10.00	40mm	931140040	\$11.50
12mm	931140012	\$7.00	22mm	931140022	\$8.50	32mm	931140032	\$10.00	45mm	931140045	\$11.50
14mm	931140014	\$7.00	24mm	931140024	\$8.50	34mm	931140034	\$10.00	50mm	931140050	\$13.00
16mm	931140016	\$7.00	26mm	931140026	\$8.50	36mm	931140036	\$10.00	55mm	931140055	\$13.00
18mm	931140018	\$8.50	28mm	931140028	\$8.50	38mm	931140038	\$10.00	60mm	931140060	\$14.00


6.5mm




6.5mm  
Fully Threaded Cancellous Screw




Standard  
Thread



3.5mm Hex Screwdriver  
(See Page 77/78)



3.2mm Pilot Hole  
(See Page 76)



Size	Code	Price	Size	Code	Price	Size	Code	Price	Size	Code	Price
30mm	931465030	\$14.00	45mm	931465045	\$16.00	55mm	931465055	\$18.00	65mm	931465065	\$18.00
35mm	931465035	\$14.00	50mm	931465050	\$17.00	60mm	931465060	\$18.00	70mm	931465070	\$19.00
40mm	931465040	\$15.00									



## Order by Phone

+44 (0) 845 045 0259  
 +1 772-562-6044  
 +61 (0) 8 9590 8850

+49 (0) 4191 8030013  
 +34 902 47 47 01  
 +39 (0) 391 11 60 784

+27 (0) 83 227 8181  
 +886 4 2258 6668

## Order Online

All of our products are available to buy online at:

[www.orthomed.co.uk/nz](http://www.orthomed.co.uk/nz)

## Screw Racks



All Orthomed Screws  
(See Page 58)

Desc	Code	Price
Small Animal Orthopaedic Tray with Screw Rack	900010000	\$690.00
As above, including full screw set	900010000SET	\$4,040.00



## Straight Plates

1.5mm

Size	Code	Price
3 hole	934104003	\$26.00
4 hole	934104004	\$28.00

Size	Code	Price
5 hole	934104005	\$30.00
6 hole	934104006	\$32.00

Size	Code	Price
20 hole (Cutttable)	935915020	\$90.00



2.0mm

Size	Code	Price
3 hole	934105003	\$23.00
4 hole	934105004	\$26.00
5 hole	934105005	\$30.00

Size	Code	Price
6 hole	934105006	\$34.00
7 hole	934105007	\$38.00
8 hole	934105008	\$42.00

Size	Code	Price
20 hole (Cutttable)	935920020	\$90.00



## Limited Contact Plates

2.0mm



Size	Code	Price
4 hole	934206004	\$32.00
5 hole	934206005	\$36.00
6 hole	934206006	\$42.00




Size	Code	Price
7 hole	934206007	\$46.00
8 hole	934206008	\$50.00



Size	Code	Price
20 hole (Cutttable)	934206020	\$108.00



**Order by Phone**

 +44 (0) 845 045 0259  
 +1 772-562-6044  
 +61 (0) 8 9590 8850

 +49 (0) 4191 8030013  
 +34 902 47 47 01  
 +39 (0) 391 11 60 784

 +27 (0) 83 227 8181  
 +866 4 2258 6668

**Order Online**

All of our products are available to buy online at:

[www.orthomed.co.uk/nz](http://www.orthomed.co.uk/nz)

67

2.7mm



Size	Code	Price
2 hole	934208002	\$52.00
3 hole	934208003	\$56.00
4 hole	934208004	\$58.00

Size	Code	Price
5 hole	934208005	\$68.00
6 hole	934208006	\$72.00
7 hole	934208007	\$76.00

Size	Code	Price
8 hole	934208008	\$82.00
9 hole	934208009	\$92.00
10 hole	934208010	\$96.00

Size	Code	Price
11 hole	934208011	\$108.00
12 hole	934208012	\$108.00

3.5mm



Size	Code	Price
4 hole	934210004	\$58.00
5 hole	934210005	\$60.00
6 hole	934210006	\$62.00

Size	Code	Price
7 hole	934210007	\$68.00
8 hole	934210008	\$72.00
9 hole	934210009	\$76.00

Size	Code	Price
10 hole	934210010	\$82.00
11 hole	934210011	\$86.00
12 hole	934210012	\$96.00

# Self Compression Plates

2.7mm



Size	Code	Price
2 hole	934108002	\$44.00
4 hole	934108004	\$48.00
5 hole	934108005	\$56.00

Size	Code	Price
6 hole	934108006	\$60.00
7 hole	934108007	\$64.00
8 hole	934108008	\$68.00

Size	Code	Price
9 hole	934108009	\$76.00
10 hole	934108010	\$80.00
11 hole	934108011	\$90.00

Size	Code	Price
12 hole	934108012	\$90.00
14 hole	934108014	\$100.00

3.5mm



Size	Code	Price
2 hole	934110002	\$44.00
3 hole	934110003	\$46.00
4 hole	934110004	\$48.00
5 hole	934110005	\$50.00

Size	Code	Price
6 hole	934110006	\$52.00
7 hole	934110007	\$56.00
8 hole	934110008	\$60.00
9 hole	934110009	\$64.00

Size	Code	Price
10 hole	934110010	\$68.00
11 hole	934110011	\$72.00
12 hole	934110012	\$80.00
13 hole	934110013	\$90.00

Size	Code	Price
14 hole	934110014	\$100.00
16 hole	934110016	\$110.00

3.5mm - Broad



Size	Code	Price
7 hole	134070000	\$100.00
8 hole	134080000	\$100.00
9 hole	134090000	\$110.00

Size	Code	Price
10 hole	134100000	\$120.00
11 hole	134110000	\$130.00
12 hole	134120000	\$140.00

Size	Code	Price
13 hole	134130000	\$160.00
14 hole	134140000	\$160.00
16 hole	134150000	\$180.00

Size	Code	Price
18 hole	134160000	\$200.00
20 hole	134170000	\$240.00
22 hole	134180000	\$280.00

68

## Order by Phone

+44 (0) 845 045 0259  
 +1 772-562-6044  
 +61 (0) 8 9590 8850

+49 (0) 4191 8030013  
 +34 902 47 47 01  
 +39 (0) 391 11 60 784

+27 (0) 83 227 8181  
 +886 4 2258 6668

## Order Online

All of our products are available to buy online at:

[www.orthomed.co.uk/nz](http://www.orthomed.co.uk/nz)

4.5<sub>mm</sub> - Narrow

Size	Code	Price
2 hole	934112002	\$56.00
3 hole	934112003	\$58.00
4 hole	934112004	\$60.00
5 hole	934112005	\$62.00
6 hole	934112006	\$64.00

Size	Code	Price
7 hole	934112007	\$66.00
8 hole	934112008	\$70.00
9 hole	934112009	\$76.00
10 hole	934112010	\$82.00
11 hole	934112011	\$90.00

Size	Code	Price
12 hole	934112012	\$95.00
13 hole	934112013	\$100.00
14 hole	934112014	\$110.00
15 hole	934112015	\$116.00
16 hole	934112016	\$130.00

Size	Code	Price
18 hole	934112018	\$136.00
20 hole	934112020	\$146.00
22 hole	934112022	\$160.00
24 hole	934112024	\$174.00

4.5<sub>mm</sub> - Broad

Size	Code	Price
6 hole	934116006	\$80.00
7 hole	934116007	\$90.00
8 hole	934116008	\$95.00
9 hole	934116009	\$100.00

Size	Code	Price
10 hole	934116010	\$105.00
11 hole	934116011	\$110.00
12 hole	934116012	\$110.00
14 hole	934116014	\$120.00

Size	Code	Price
16 hole	934116016	\$130.00
18 hole	934116018	\$140.00
20 hole	934116020	\$150.00
22 hole	934116022	\$160.00

Size	Code	Price
24 hole	934116024	\$170.00
26 hole	934116026	\$180.00

## T Plates

Size	Code	Price
2.0mm	935920000	\$20.00
2.7mm	935918000	\$20.00

3.5<sub>mm</sub> - Right Angled









Size	Code	Price
3HHX3HT	935913033	\$50.00
3HHX4HT	935913034	\$50.00
3HHX5HT	935913035	\$60.00

Size	Code	Price
4HHX4HT	935913044	\$80.00
4HHX6HT	935913046	\$80.00





**Order by Phone**

 +44 (0) 845 045 0259   
  +49 (0) 4191 8030013   
  +27 (0) 83 227 8181  
 +1 772-562-6044   
  +34 902 47 47 01   
  +886 4 2258 6668  
 +61 (0) 8 9590 8850   
  +39 (0) 391 11 60 784

**Order Online**

All of our products are available to buy online at:

[www.orthomed.co.uk/nz](http://www.orthomed.co.uk/nz)



**3.5mm - Oblique Angled**

Size	Code	Price
<b>3HHX3HT</b>	935913133	<b>\$54.00</b>
<b>3HHX4HT</b>	935913134	<b>\$54.00</b>

Size	Code	Price
<b>3HHX5HT</b>	935913135	<b>\$54.00</b>



**3.5mm - Thick Oblique Angled**

Size	Code	Price
<b>3HHX3HT</b>	935914133	<b>\$68.00</b>
<b>3HHX4HT</b>	935914134	<b>\$68.00</b>

Size	Code	Price
<b>3HHX5HT</b>	935914135	<b>\$75.00</b>



# Acetabular Plates

Size	Code	Price
<b>2.0mm</b>	107401151	<b>\$80.00</b>
<b>2.7mm</b>	107401150	<b>\$80.00</b>

Size	Code	Price
<b>3.5mm</b>	107401149	<b>\$80.00</b>



# Cortical Nuts

Size	Code	Price
<b>3.5mm</b>	932999010	<b>\$14.00</b>
<b>4.5mm</b>	932999011	<b>\$14.00</b>



# Washers




**Standard**




Size	Code	Price
<b>1.5-2.0mm</b>	932999003	<b>\$5.00</b>
<b>2.7-4.0mm</b>	932999002	<b>\$5.00</b>
<b>4.5-6.5mm</b>	932999001	<b>\$8.00</b>





70

## Order by Phone

 +44 (0) 845 045 0259  
 +1 772-562-6044  
 +61 (0) 8 9590 8850

 +49 (0) 4191 8030013  
 +34 902 47 47 01  
 +39 (0) 391 11 60 784

 +27 (0) 83 227 8181  
 +886 4 2258 6668

## Order Online

All of our products are available to buy online at:

[www.orthomed.co.uk/nz](http://www.orthomed.co.uk/nz)

## Ligament

Size	Code	Price
2.7mm	932999014	\$40.00
3.5mm	932999013	\$52.00
4.5-6.0mm	932999012	\$56.00



## Cerclage Wire

10<sub>m</sub> - Coil

Size	Code	Price
0.6mm	938506010	\$26.00
0.8mm	938508010	\$36.00
1.0mm	938510010	\$46.00
1.2mm	938512010	\$56.00
1.5mm	938515010	\$80.00

280<sub>mm</sub> - Loop

Size	Code	Price
0.8mm	938408280	\$8.00
1.0mm	938410280	\$8.00



Size	Code	Price
1.2mm	938412280	\$8.00
1.5mm	938415280	\$8.00

## Steinmann Pins

300<sub>mm</sub>

Single Pointed Tip









Size	Code	Price
1.6mm	938016300	\$8.00
1.8mm	938018300	\$8.00
2.0mm	938020300	\$8.50
2.4mm	938024300	\$10.00



Size	Code	Price
2.8mm	938028300	\$11.00
3.2mm	938032300	\$12.00
3.6mm	938036300	\$13.50
4.0mm	938040300	\$17.60

Size	Code	Price
4.8mm	938048300	\$18.40
6.0mm	938060300	\$26.50

## Order by Phone

 +44 (0) 845 045 0259  
 +1 772-562-6044  
 +61 (0) 8 9590 8850

 +49 (0) 4191 8030013  
 +34 902 47 47 01  
 +39 (0) 391 11 60 784

 +27 (0) 83 227 8181  
 +886 4 2258 6668

## Order Online

All of our products are available to buy online at:

[www.orthomed.co.uk/nz](http://www.orthomed.co.uk/nz)

71

# Trochar Pins (Pack of 10)

## 150<sub>mm</sub> - Single Point



Single Pointed Tip



Size	Code	Price
0.8mm	937508150	\$32.00
1.0mm	937510150	\$32.00

Size	Code	Price
1.2mm	937512150	\$34.00
1.6mm	937516150	\$34.00

Size	Code	Price
2.0mm	937520150	\$34.00

## 150<sub>mm</sub> - Double Point



Double Pointed Tip



Size	Code	Price
0.8mm	937608150	\$42.00
1.0mm	937610150	\$42.00

Size	Code	Price
1.2mm	937612150	\$45.00
1.6mm	937616150	\$45.00

Size	Code	Price
2.0mm	937620150	\$45.00







# Instruments

---

## 07

Our reputation for premium quality is ably demonstrated by our range of core orthopaedic instrumentation, with special mention due to the bespoke instruments that is a fundamental part of orthopaedic procedures. All our instruments are manufactured with the same attention to detail found across the rest of our range, from premium quality materials, all adhering closely to current manufacturing standards.




### In This Section




Metzenbaum Scissors	74	Taps	76
Olsen Needle Holders	74	Counter Sinks	77
Spencer Wells Forceps	74	Quick Coupling Handles	77
Mayo Scissors	74	HEX Screwdrivers	77
Allis Tissue Forceps	75	Depth Gauges	78
Reduction Forceps	75	Drill Guides	78
Dissecting Forceps	75	Drill Sleeve	79
Shardel Towel Clip	75	Bending Irons	79
Wire Cutters	76	Retractors	79
Drill Bits	76	Periosteal Elevators	79



---

74

## Order by Phone

 +44 (0) 845 045 0259  
 +1 772-562-6044  
 +61 (0) 8 9590 8850

 +49 (0) 4191 8030013  
 +34 902 47 47 01  
 +39 (0) 391 11 60 784

 +27 (0) 83 227 8181  
 +886 4 2258 6668

## Order Online

All of our products are available to buy online at:

[www.orthomed.co.uk/nz](http://www.orthomed.co.uk/nz)

## Metzenbaum Scissors

Size	Code	Price
145mm	910001631	\$105.00
180mm	910001632	\$125.00

Size	Code	Price
180mm Tungsten Carbide	910002401TC	\$170.00



## Olsen Needle Holders

Size	Code	Price
145mm	910002370	\$138.00

Size	Code	Price
145mm Tungsten Carbide	910002371TC	\$210.00
170mm Tungsten Carbide	910002400TC	\$217.00



## Spencer Wells Forceps

Size	Code	Price
125mm	910001593	\$77.00


Size	Code	Price
180mm	910001595	\$90.00




## Mayo Scissors


Size	Code	Price
140mm	910001627	\$92.50

Size	Code	Price
165mm	910001439	\$108.00




## Order by Phone

 +44 (0) 845 045 0259

 +1 772-562-6044


 +61 (0) 8 9590 8850

 +49 (0) 4191 8030013

 +34 902 47 47 01

 +39 (0) 391 11 60 784

 +27 (0) 83 227 8181

 +886 4 2258 6668

## Order Online

All of our products are available to buy online at:

[www.orthomed.co.uk/nz](http://www.orthomed.co.uk/nz)

75

## Allis Tissue Forceps

Size	Code	Price
2/3 teeth	910001699	\$85.00
3/4 teeth	910001700	\$85.00

Size	Code	Price
4/5 teeth	910001701	\$114.00



## Reduction Forceps

Size	Code	Price
Small	900500004	\$310.00

Size	Code	Price
Large	900500008	\$325.00



## Dissecting Forceps

Gillies  
Non-toothedMcIndoe  
Non-toothedWaugh's  
1 in 2 TeethWaugh's  
Non-toothed

Type	Code	Price
Gillies Non-toothed	910001651	\$62.50
McIndoe Non-toothed	910001658	\$49.50

Type	Code	Price
Waugh's 1 into 2 teeth	910001659	\$55.50
Waugh's Non-toothed	910001661	\$60.00






## Shardel Towel Clip




Size	Code	Price
90mm	910002372	\$35.50





76

## Order by Phone

 +44 (0) 845 045 0259  
 +1 772-562-6044  
 +61 (0) 8 9590 8850

 +49 (0) 4191 8030013  
 +34 902 47 47 01  
 +39 (0) 391 11 60 784


 +27 (0) 83 227 8181  
 +886 4 2258 6668

## Order Online

All of our products are available to buy online at:

[www.orthomed.co.uk/nz](http://www.orthomed.co.uk/nz)

## Wire Cutters



Size	Code	Price
180mm	900414002	\$410.00

Size	Code	Price
215mm	900414004	\$620.00



## Drill Bits

### Jacobs Chuck



Jacobs Chuck



Use Drill  
(See Page 82)



Size	Code	Price
1.1mm	900000011	\$40.00
1.5mm	900000015	\$40.00
1.8mm	900000018	\$40.00

Size	Code	Price
2.0mm	900000020	\$40.00
2.5mm	900000025	\$50.00
2.7mm	900000027	\$50.00

Size	Code	Price
3.2mm	900000032	\$50.00
3.5mm	900000035	\$60.00
4.0mm	900000040	\$60.00

Size	Code	Price
4.5mm	900000045	\$60.00

### Quick Coupling



Quick Coupling



Use Quick Coupling Handles  
(See Page 77)



Size	Code	Price
1.1mm	900100011	\$44.00
1.5mm	900100015	\$44.00

Size	Code	Price
2.0mm	900100020	\$44.00
2.5mm	900100025	\$55.00

Size	Code	Price
2.7mm	900100027	\$55.00
3.2mm	900100032	\$55.00

Size	Code	Price
3.5mm	900100035	\$65.00
4.5mm	900100045	\$65.00

## Taps

### Cortical



Quick Coupling



Use Quick Coupling Handles  
(See Page 77)



Use Cortical Screws  
(See Page 58)






Size	Code	Price
1.5mm	900200150	\$116.00
2.0mm	900200200	\$116.00




Size	Code	Price
2.7mm	900200270	\$116.00
3.5mm	900200350	\$116.00



Size	Code	Price
4.5mm	900200450	\$116.00



## Order by Phone

 +44 (0) 845 045 0259  
 +1 772-562-6044  
 +61 (0) 8 9590 8850

 +49 (0) 4191 8030013  
 +34 902 47 47 01  
 +39 (0) 391 11 60 784

 +27 (0) 83 227 8181  
 +886 4 2258 6668

## Order Online

All of our products are available to buy online at:

[www.orthomed.co.uk/nz](http://www.orthomed.co.uk/nz)

77

## Cancellous



Quick Coupling

Use Quick Coupling Handles  
(See Page 77)Use Cancellous Screws  
(See Page 65)

Size	Code	Price
4.0mm	900201400	\$116.00

Size	Code	Price
6.5mm	900201650	\$116.00

## Counter Sinks



Quick Coupling

Use Quick Coupling Handles  
(See Page 77)

Size	Code	Price
1.5-2.0mm	900415001	\$105.00
2.7mm	900427001	\$105.00

Size	Code	Price
3.5mm	900435001	\$105.00
4.5mm (with handle)	900445001	\$168.00



## Quick Coupling Handles

For Quick Coupling Instruments  
(See Page 76/77)

Type	Code	Price
Mini Straight	900420001	POA
Pear Shape	900420002	POA

Type	Code	Price
T Shape	900420003	\$199.00



## HEX Screwdrivers

## Quick Coupling Shafts

Hexagonal  
HeadUse Quick Coupling Handles  
(See Page 77)




Size	Code	Price
1.5mm	900425005 (For 1.5-2.4mm Screws)	\$115.00




Size	Code	Price
2.5mm	900425030 (For 2.7-4.0mm Screws)	\$115.00



Size	Code	Price
3.5mm	900425050 (For 4.5-6.5mm Screws)	\$115.00

78

**Order by Phone**

 +44 (0) 845 045 0259  
 +1 772-562-6044  
 +61 (0) 8 9590 8850

 +49 (0) 4191 8030013  
 +34 902 47 47 01  
 +39 (0) 391 11 60 784

 +27 (0) 83 227 8181  
 +886 4 2258 6668

**Order Online**

All of our products are available to buy online at:

[www.orthomed.co.uk/nz](http://www.orthomed.co.uk/nz)

**Complete Soft Grip**



Hexagonal Head



Size	Code	Price
<b>1.5mm</b>	900425060 (For 1.5-2.4mm Screws)	<b>\$336.00</b>

Size	Code	Price
<b>2.5mm</b>	900425065 (For 2.7-4.0mm Screws)	<b>\$336.00</b>

Size	Code	Price
<b>3.5mm</b>	900425070 (For 4.5-6.5mm Screws)	<b>\$336.00</b>

## Depth Gauges

Size	Code	Price
<b>1.5/2.0mm</b>	900400152	<b>\$600.00</b>
<b>2.7/4.0mm</b>	900400274	<b>\$600.00</b>



## Drill Guides

**Double**




For Drill Bits,  
(See Page 76)

Size	Code	Price
<b>1.1/1.5mm</b>	900200115	<b>\$250.00</b>
<b>1.5/2.0mm</b>	900200152	<b>\$250.00</b>
<b>2.0/2.7mm</b>	900200207	<b>\$250.00</b>

Size	Code	Price
<b>2.5/3.5mm</b>	900200255	<b>\$250.00</b>
<b>3.2/4.5mm</b>	900200324	<b>\$250.00</b>



**Compression**




For Drill Bits,  
(See Page 76)


Size	Code	Price
<b>2.0mm</b>	900320000	<b>\$273.00</b>
<b>2.7mm</b>	900327000	<b>\$273.00</b>

Size	Code	Price
<b>3.5mm</b>	900335000	<b>\$283.50</b>
<b>4.5mm</b>	900345000	<b>\$283.50</b>





## Order by Phone


 +44 (0) 845 045 0259


 +1 772-562-6044


 +61 (0) 8 9590 8850

 +49 (0) 4191 8030013

 +34 902 47 47 01

 +39 (0) 391 11 60 784

 +27 (0) 83 227 8181

 +886 4 2258 6668

## Order Online

All of our products are available to buy online at:

[www.orthomed.co.uk/nz](http://www.orthomed.co.uk/nz)

79

## Drill Sleeve

Size	Code	Price
3.5/2.5mm	900303525	\$84.00



## Bending Irons

Size	Code	Price
1.5-2.0mm	900411015	\$147.00
2.7-4.0mm (Left)	900411127	\$157.50

Size	Code	Price
2.7-4.0mm (Right)	900411227	\$157.50



## Retractors



6mm

8mm

15mm

Size	Code	Price
6mm	900530006	\$105.00
8mm	900530008	\$105.00

Size	Code	Price
15mm	900530015	\$147.00



## Periosteal Elevators

Size	Code	Price
3mm	900535003	\$205.00

Size	Code	Price
6mm	900535006	\$215.00







# Power Tools

## Kits & Spares


---

08


Our high performance, premium quality saw and drill come supplied as a complete kit, comprising of 2 batteries, charger, sterile shroud and chuck.

The saw also comes complete with blades – ideal for use with our MMP system. Both come with a two year guarantee.

### In This Section






Drill








Saw

---

## Order by Phone

 +44 (0) 845 045 0259  
 +1 772-562-6044  
 +61 (0) 8 9590 8850

 +49 (0) 4191 8030013  
 +34 902 47 47 01  
 +39 (0) 391 11 60 784

 +27 (0) 83 227 8181  
 +886 4 2258 6668


## Order Online

All of our products are available to buy online at:

[www.orthomed.co.uk/nz](http://www.orthomed.co.uk/nz)

# Bone Drill Set




Size	Code	Price
Bone Drill Set	OK8002	\$2,800.00








Unit rise in temperature	≤25°C
Charge	100V-250V; 50/60Hz
Battery Voltage	7.2V

Charging Time	About 2 hours
Non-load noise	≤40cB(A)
Can be sterilized up to	135°C

**Order by Phone**

 +44 (0) 845 045 0259  
 +1 772-562-6044  
 +61 (0) 8 9590 8850

 +49 (0) 4191 8030013  
 +34 902 47 47 01  
 +39 (0) 391 11 60 784

 +27 (0) 83 227 8181  
 +886 4 2258 6668

**Order Online**

All of our products are available to buy online at:

[www.orthomed.co.uk/nz](http://www.orthomed.co.uk/nz)



# Sagittal Saw Set

Size	Code	Price
Sagittal Saw Set	OK8001	\$2,800.00

Blade Locking Position  
 Blade Oscillation Arc 5°  
 Operating Speed: 0~18000 rpm

**Sagittal Saw Blades**



**Sagittal Saw Set**



# Sagittal Saw Blades

Size	Code	Price
MMP Saw Blade 40 x 12 x 0.6mm	OKMMP001	\$28.00
Saw Blade 40 x 8 x 0.6mm	OKMMP002	\$28.00
Saw Blade 40 x 5 x 0.6mm	OKMMP003	\$28.00



Unit rise in temperature	≤25°C
Charge	100V-250V; 50/60Hz
Battery Voltage	7.2V

Charging Time	About 2 hours
Non-load noise	≤70cB(A)
Features can be sterilized up to	135°C







# More From Orthomed

---



We're passionate about the industry and the people that work in it, so we want to be able to offer the best possible support to all of our customers.

We've come up with a range of easily accessible tutorial material to offer advice before, during and after surgery.

We're a friendly bunch at Orthomed, so if you want to chat rather than just read or watch, we're only at the other end of the phone and we'd be happy to help.



Order by Phone

+44 (0) 845 045 0259  
+1 772-562-6044  
+61 (0) 8 9590 8850

+49 (0) 4191 8030013  
+34 902 47 47 01  
+39 (0) 391 11 60 784

+27 (0) 83 227 8181  
+886 4 2258 6668

Order Online

All of our products are available to buy online at:

[www.orthomed.co.uk/nz](http://www.orthomed.co.uk/nz)

# Supporting your journey in orthopaedics

## Expertise you can trust at every step

Orthomed are not just about selling orthopaedic products, we are highly committed to helping develop the industry and the people that work and excel in it.



**Order by Phone**

🇬🇧 +44 (0) 845 045 0259  
 🇺🇸 +1 772-562-6044  
 🇦🇺 +61 (0) 8 9590 8850

🇩🇪 +49 (0) 4191 8030013  
 🇪🇸 +34 902 47 47 01  
 🇮🇹 +39 (0) 391 11 60 784

🇳🇿 +27 (0) 83 227 8181  
 🇵🇰 +886 4 2258 6668

**Order Online**

All of our products are available to buy online at:

[www.orthomed.co.uk/nz](http://www.orthomed.co.uk/nz)

87

# ourney



We take education very seriously and that's why we strive to offer the highest quality hands-on practical workshops, training materials, & on-going support to all our customers.




**Customers benefit from:**




- ✔ Dry Lab Workshops
- ✔ Wet Lab Workshops
- ✔ Technical Support
- ✔ Webinars
- ✔ Userguides
- ✔ Templates
- ✔ Case Studies
- ✔ Videos
- ✔ Animations
- ✔ Expert Advice
- ✔ Catalogues
- ✔ Surgical Support



Visit our website for more info [www.orthomed.co.uk/nz](http://www.orthomed.co.uk/nz)



#### Order by Phone

 +44 (0) 845 045 0259  
 +1 772-562-6044  
 +61 (0) 8 9590 8850

 +49 (0) 4191 8030013  
 +34 902 47 47 01  
 +39 (0) 391 11 60 784

 +27 (0) 83 227 8181  
 +886 4 2258 6668

#### Order Online

All of our products are available to buy online at:

[www.orthomed.co.uk/nz](http://www.orthomed.co.uk/nz)

# Orthomed around the world

## Supporting all of our customers

Although Orthomed was born and bred in the UK (and is very proud of its Britishness!), our desire to bring innovative implants to the wider veterinary community has resulted in swift expansion across the globe.

Orthomed currently has sister companies in the USA, Australia, Germany, Italy, Spain and Portugal, and South Africa, and this is just the start of our global expansion plans.

Overseas customers do not need to visit the UK to take advantage of our educational offerings – they will find themselves just as well catered for as their British counterparts, with workshops taking place in each of these countries. Although, they will of course be warmly welcomed if they choose to visit the UK and take part in a workshop on Orthomed's home soil...



### Orthomed (UK) Ltd

Majestic House  
29 Green Street, Huddersfield  
West Yorkshire, HD1 5DQ  
United Kingdom

**Tel:** +44 (0) 845 045 0259  
**Fax:** +44 (0) 845 603 2456  
**Email:** [info@orthomed.co.uk](mailto:info@orthomed.co.uk)



### Orthomed North America Inc.

1179 18th Place  
Vero Beach  
Florida, 32960  
USA

**Tel:** +1 772-562-6044  
**Fax:** +1 772-562-6046  
**Email:** [info@orthomedusa.com](mailto:info@orthomedusa.com)



### Orthomed Technology GmbH

Am Schaafredder 17  
24568, Kaltenkirchen  
Germany

**Tel:** +49 (0) 4191 8030013  
**Fax:** +49 (0) 4191 8030014  
**Email:** [info@orthomedeu.com](mailto:info@orthomedeu.com)






### Orthomed (SA) Pty Ltd




90 Henry Rd  
Shere  
Pretoria, 0084  
South Africa



**Tel:** +27 (0) 83 227 8181  
**Fax:** +27 (0) 86 683 2511  
**Email:** [info@orthomedsa.co.za](mailto:info@orthomedsa.co.za)



**Order by Phone**

 +44 (0) 845 045 0259  
 +1 772-562-6044  
 +61 (0) 8 9590 8850

 +49 (0) 4191 8030013  
 +34 902 47 47 01  
 +39 (0) 391 11 60 784

 +27 (0) 83 227 8181  
 +886 4 2258 6668

**Order Online**

All of our products are available to buy online at:

[www.orthomed.co.uk/nz](http://www.orthomed.co.uk/nz)



[www.orthomed.co.uk/nz](http://www.orthomed.co.uk/nz)



**Orthomed Australasia Pty Ltd**

PO Box 881  
 Mandurah  
 WA, 6210  
 Australia

**Tel:** +61 8 9590 8850  
**Fax:** +61 8 9510 9001  
**Email:** [info@orthomed.com.au](mailto:info@orthomed.com.au)



**B Braun VetCare S.A.**

*Authorised distributor for Spain and Portugal*

Carretera de Terrassa, 121  
 08191, Rubi  
 Barcelona  
 Spain

**Tel:** +34 902 47 47 01  
**Fax:** +34 902 48 48 01  
**Email:** [vet.es@bbraun.com](mailto:vet.es@bbraun.com)



**Orthomed IT**

Via Crociata 18  
 35028, Piove di Sacco  
 Padova  
 Italy

**Tel:** +39 (0) 391 11 60 784  
**Email:** [barbara@orthomed-italia.it](mailto:barbara@orthomed-italia.it)



**Orthomed Taiwan, EPOCH Ltd**

No.41, Sec. 2, Huizhong Rd  
 Xitun Dist.  
 Taichung City 40756  
 Taiwan

**Tel:** +886 4 2258 6668  
**Fax:** +886 4 2258 6558  
**Email:** [info@ahoss.com.tw](mailto:info@ahoss.com.tw)



## Orthomed (UK) Ltd

Majestic House  
29 Green Street, Huddersfield  
West Yorkshire, HD1 5DQ  
United Kingdom

**Tel:** +44 (0) 845 045 0259

**Fax:** +44 (0) 845 603 2456

**Email:** info@orthomed.co.uk



## Orthomed (SA) Pty Ltd

90 Henry Rd  
Shere  
Pretoria, 0084  
South Africa

**Tel:** +27 (0) 83 227 8181

**Fax:** +27 (0) 86 683 2511

**Email:** info@orthomedsa.co.za



## Orthomed North America Inc.

1179 18th Place  
Vero Beach  
Florida, 32960  
USA

**Tel:** +1 772-562-6044

**Fax:** +1 772-562-6046

**Email:** info@orthomedusa.com



## Orthomed Australasia Pty Ltd

PO Box 881  
Mandurah  
WA, 6210  
Australia

**Tel:** +61 (0) 8 9590 8850

**Fax:** +61 (0) 8 9510 9001

**Email:** info@orthomed.com.au



## Orthomed Technology GmbH

Am Schaafredder 17  
24568, Kaltenkirchen  
Germany

**Tel:** +49 (0) 4191 8030013

**Fax:** +49 (0) 4191 8030014

**Email:** info@orthomedeu.com



## Orthomed IT

Via Crociata 18  
35028, Piove di Sacco  
Padova  
Italy

**Tel:** +39 (0) 391 11 60 784

**Email:** barbara@orthomed-italia.it



*Authorised distributor for Spain and Portugal*

## B Braun VetCare S.A.

Carretera de Terrassa, 121  
08191, Rubi  
Barcelona  
Spain

**Tel:** +34 902 47 47 01

**Fax:** +34 902 48 48 01

**Email:** vet.es@bbraun.com



## Orthomed Taiwan, EPOCH Ltd

No.41, Sec. 2, Huizhong Rd.,  
Xitun Dist.  
Taichung City 40756  
Taiwan

**Tel:** +886 4 2258 6668

**Fax:** +886 4 2258 6558

**Email:** info@ahoss.com.tw



The home of  
**Veterinary Orthopaedics**

[www.orthomed.co.uk/nz](http://www.orthomed.co.uk/nz)